

Frozen shipping made simple

Need to safely and cost-effectively ship your high-end, temperature-sensitive biogenetics and biomaterials to destinations around the globe? Here's an easy, worry-free solution.

The FedEx® Deep Frozen Shipping **Solution** couples innovative liquid nitrogen dry vapor technology with end-to-end shipment monitoring and visibility for superior integrity and security for your high-value frozen shipments worldwide.

Peace of mind

The FedEx Deep Frozen Shipping Solution offers you several advantages over dry-ice shipping, including:

- Longer hold times. The solution provides you with up to 10 days of holding time (below -150°C) to support complex international moves and eliminate sublimation and the need for re-icing that's required with dry ice shipments.
- Regulatory compliance. The non-hazardous dry shipper meets International Air Transport Association (IATA) requirements for both infectious substances (UN3373) and non-infectious clinical samples. You can also ship biomaterials that are classified as dangerous goods.
- Inventory elimination. You won't need to purchase, inventory or manage insulated boxes and dry ice or train personnel on handling hazardous materials.
- Advanced tracking, monitoring and care. With every shipment, you automatically receive the FedEx Priority Alert™ service that provides proactive monitoring and intervention, if needed, from origin to destination.

With the FedEx® Deep Frozen Shipping Solution, you get the most sophisticated technology for your frozen shipments in just a few easy steps:

Place your order customized, secure FedEx portal.



FedEx picks up the empty container and returns it to a recharging station.



The recipient unpacks the materials and prepares return shipment of empty container using the label provided.



FedEx delivers the dry online through a · · · · shipping container, pre-charged, to your pickup point.



Load your materials into the container and apply the shipping label provided.



FedEx picks up the shipment and delivers it to the recipient.



Safe, reliable and earth-friendly

Enjoy all the other advantages that the innovative FedEx Deep Frozen Shipping Solution has to offer, including:

- Reliability. Receive proactive shipment monitoring from FedEx every step of the way to ensure end-toend security and integrity.
- Flexibility. Choose the FedEx Express[®] service that best suits your domestic delivery needs from urgent to less critical.
- Visibility. Enter your order and track the shipment anywhere in the world via a customized, secure FedEx portal.
- Earth-friendly. Reduce impact on the environment as the dry shipper is recycled with negligible and fully recyclable cardboard waste and liquid nitrogen evaporates harmlessly.

	Traditional dry-ice shipments	FedEx® Deep Frozen Shipping Solution
Container considered non-dangerous goods and eliminates need for licensed handler?	No	✓
Includes pre-arranged shipment logistics?	No	✓
Has up to 10 days of holding time to accommodate international travel?	No	✓
Requires no inventory to buy and maintain?	No	✓
Exempted from quantity restrictions on aircraft?	No	✓
Made of recyclable materials?	No	✓

Variety of sizes

Choose the dry shipper that best suits your needs:

- **Standard dry shipper,** for quantities of approximately 75 1.8 or 2.0 milliliter vials.
- Large-volume dry shipper, for quantities of approximately 500 1.8 or 2.0 milliliter vials.

	Standard dry shipper	Large-volume dry shipper
Specimen chamber		
Diameter	2.8"	8.5"
Depth	11.5"	12"
Storage capacity	75 of 1.8 and 2.0 ml vials	5 cryovial boxes (5"x 5"x 2") –
	15 canes with 5 vial clips	81 2.0 ml vials per box
	75 cryovials, 2.0 ml	5 cryovial boxes (5.25"x 5.25"x 2") – 100 2.0 ml vials per box
		6 Tyvek bags – 15 5-vial canes per bag
Weight		
Dry	15 lbs.	50 lbs.
Wet	20 lbs.	71 lbs.
Domestic dimensional	21.07 lbs.	79.68 lbs.
International dimensional	24.87 lbs.	94.01 lbs.
Outer box	12" x 24"	22" x 27"



More info

- Contact your FedEx account executive
- Go to fedex.com/healthcare
- Call 1.877.726.2091